

MULTI FUNCTIONAL HAND-PIECES

LARGE	Abdomen and Back
MEDIUM	Arm and Thigh
SMALL	Face and Neck

SYSTEM SPECIFICATIONS

Radio Frequency	0.8 MHz
Power	70 W
Vacuum	200 mmHg
LED Wavelength	640/460 nm
Rating	100-240 V, 50/60 Hz
Weight	7.2 kg
Dimension	400 (D) x 405 (W) x 220 (H) mm



All In One Solution
by Cluederm.
Expose Your Potential.

RTB012016

-
- Body Shaping
- Cellulite Reduction
- Facial Lifting & Tightening
- Wrinkle Reduction
-



RIDM
Technology



SAS
Technology



Cluederm™ LIFTING PATCH | BODY PATCH



Entire Beauty Solution
Cluederm™

info@cluederm.com | www.cluederm.com

Images and texts are intellectual property of Cluederm. Copying of this material can be subject to charges of both civil and criminal law of legal justice.
Copyright to Cluederm © 2016

Feel the Rhythm
Refit.™

Smart Facial Rejuvenation and Body Shaping System

Cluederm™
Entire Beauty Solution

Feel the Rhythm


Refit™

Smart Facial Rejuvenation
and Body Shaping System



• RIDM (Redesign Integrated Dynamic Movement) Technology

Addressing the issues of operator fatigue and complicated procedures, treatment procedures are easier to operate with the support of RIDM technology. There is no need to alter the parameters or check the time during treatment procedures on a regular basis.



• SAS (Safe Anti-Spark) Technology


RF energy controls, SAS Technology, have been set in place to ensure the safety and protection of the patient from possible sparking. This commonly occurs with other similar devices utilizing RF technology, due to unstable power levels. The Refit equipped with SAS technology allows physicians to operate with confidence when using higher frequency levels.



User orientated. Patient orientated. Result orientated.

• Integrated Technology

The synergy of multi-polar RF, Vacuum suction, and LED technology competently achieve lifting and tightening effects, cellulite reduction and body contouring. The combination of these three functions allow users to expose their potential with minimal pain and no downtime.



Multi-polar RF

The RF energy delivers thermal energy into the dermal and subcutaneous layer, promoting faster metabolism, tightening collagen fibers and improving skin elasticity.

Vacuum Suction


The vacuum suction increases blood circulation, reduces the adhesive contact between fat cell clusters, and therefore reshapes the body. Additionally, it assists in effective penetration of RF energy delivery to deeper levels of the skin.

LED

Custom selected wavelength of 640nm and 460nm works to increase blood circulation and improve skin tone.

• Dynamic Pulsation


Optional selection of pulse levels are available to customize treatment procedures with minimal noise levels.



• Multiple Hand-pieces

Hand-pieces are available in 3 sizes (Large, Medium and Small).

- The smaller hand-piece is used to treat the face and neck.
- The medium hand-piece is used to treat the arm and thigh.
- The large hand-piece is used to treat the back and abdomen.



• Application

Abdomen



Thigh



Arm



Face

