


Feel the Rhythm

Refit.™



What is Refit treatment?

The Refit, the latest in Radio Frequency technology represents the synergy between Multi-Polar RF, vacuum suction and LED functions to competently achieve collagen contraction and production, facial lifting and tightening, cellulite reduction, wrinkle reduction, and skin rejuvenation treatments. RF technology since its introduction in early 2000, has been mastered to expose your potential with zero pain, no downtime, and entirely at your convenience.



FAQ

Does the Refit treatment hurt?

Pain VS Reward is a crucial factor when assessing the right treatment for patients. What if there was no pain? The Refit provides a solution with no pain, no downtime and no additional concerns. The adjustable vacuum suction creates a shield for any discomfort experienced by the emission of the electrical current focused on the tissues beneath the surface of the skin.

How long does the treatment take?

The treatment procedure duration is estimated to be between 10-20 minutes and varies based on your treatment plan. The RF, Vacuum and Pulse levels are adjustable based around your preference which will determine the length of your treatment time.

Who is a suitable candidate for Refit treatment?

The Refit and its technology support all skin types and achieves skin laxity, cellulite reduction, collagen development, lifting and tightening procedures, shaping your body to accomplish your true potential. As the Refit uses radio frequency technology, candidates who have implants, internal defibrillators or pacemakers are not suitable for this treatment type. For other contraindications please inquire at your consultation session.

How does the Refit work?

The Refit utilizes radio frequency energy to effectively tighten and rejuvenate the skin by manipulating the tissue's resistance within the various layers of the skin. The RF energy is transformed into thermal energy to trigger collagen development and to achieve skin rejuvenation.

Clinical Cases



Multi-polar RF

The RF energy delivers thermal energy into the dermal and subcutaneous layer, promoting faster metabolism, tightening collagen fibers and improving skin elasticity.

Vacuum Suction

The vacuum suction increases blood circulation, reduces the adhesive contact between fat cell clusters, and therefore reshapes the body. Additionally, it assists in effective penetration of RF energy delivery to deeper levels of the skin.

LED

The variation of LED provides vitality to dull skin tones, increases cell regeneration and renewal, combats wrinkles and triggers collagen production, to achieve a smooth and balanced complexion.



SafeAnti-SparkTechnology

Application

ABDOMEN



THIGH



ARM



FACE




A global leader in aesthetics
Entire Beauty Solution
by Cluederm

Body Shaping
Cellulite Reduction
Facial Lifting & Tightening
Wrinkle Reduction

For more information visit our website or email us at
www.cluederm.com info@cluederm.com

Feel the Difference



Cluederm™ LIFTING PATCH | BODY PATCH



Cluederm™
Entire Beauty Solution

Feel the Rhythm™
Refit™

Smart Facial Rejuvenation and Body Shaping System

Cluederm™
Entire Beauty Solution